

Hybrid High Energy Vibratory Ball Mill

Item Number: KT-BM400



Introduction

KT-BM400 is used for rapid grinding or mixing of dry, wet and frozen small amount of samples in the laboratory. It can be configured with two 50ml ball mill jars

[Learn More](#)

Recommended maximum injection size	< 8 mm
Sample particle size range	~ 5um
Grinding method	Dry Grinding Wet Grinding Cryogenic Grinding
Grinding platform (number of cans)	2
Movement method	plane arc reciprocating motion
Vibration frequency	180-1800 r/min
Control method	LCD display
Time setting	59:59:59 (h/m/s)
Ball mill material	Tungsten carbide zirconia stainless steel MC nylon PTFE, etc.
Ball mill tank volume	50ml 30ml
The sealing method of the ball mill jar	Sealing ring + screw fastening
Grinding ball size	1-25m optional
Grinding ball material	Tungsten carbide Zirconium oxide Agate Stainless steel
Cell breaking adapter volume	10x0.2ml 5x2ml 5x5ml
Electrical description	100-120V/200-240V AC, 50-60Hz, 150W
Power port	National Standard European Standard American Standard British Standard etc.
Net weight	48kg
Dimensions (LWH)	470*360*250
Protection class	IP63
Standard	CE