

Assemble Lab Cylindrical Press Mold

Item Number: PMAC



Introduction

Get reliable and precise molding with Assemble Lab Cylindrical Press Mold. Perfect for ultra-fine powder or delicate samples, widely used in material research and development.

[Learn More](#)

| Model | PMAC-A | PMAC-B | PMAC-C | PMAC-D | PMAC-E |
|---------------------|-----------------------|------------------------|--------------------|----------------------|------------------------|
| Material | ASSAB+17 | Cr12MoV | | | |
| Indenter hardness | HRC68-HRC70 | HRC60-HRC62 | | | |
| Sample size | φ3/φ4/φ5/φ6/φ8/φ10 mm | φ12/φ13/φ15/φ18/φ20 mm | φ30/φ40 mm | φ50/φ60 mm | φ70/φ80 mm |
| Cavity depth | 30 mm | 40 mm | 50 mm | 55 mm | 60 mm |
| External dimensions | φ43×93 mm | φ53×120 mm | φ73×133/φ95×133 mm | φ115×150/φ127×150 mm | φ153×180mm/φ180×180 mm |
| Weight(Kg) | 0.75 | 1.2 | 3.8/6.3 | 14/20 | 30/40 |

Other sizes can be customized