


Assemble Square Lab Press Mold

Item Number: PMAS



Introduction

Achieve perfect sample preparation with Assemble Square Lab Press Mold. Quick disassembly eliminates sample deformation. Perfect for battery, cement, ceramics, and more. Customizable sizes available.

[Learn More](#)

Model	PMAS-A	PMAS-B	PMAS-C	PMAS-D	PMAS-E
Material	Cr12MoV				
Indenter hardness	HRC60-HRC62				
Sample size	3×3 /4×4 /5×5 /6×6 /8×8 /10×10 mm	12×12 /15×15 /18×18 /20×20 mm	30×30 /40×40 mm	50×50 /60×60 mm	70×70 /80×80 mm
Cavity depth	30 mm	40 mm	50 mm	55 mm	60 mm
External dimensions	φ53×120 mm	φ73×133 mm	φ95×133/φ115×133 mm	φ127×150/φ153×150 mm	φ180×180mm/φ200×180 mm
Weight(Kg)	1.2	3.6	7/14	20/30	40/50

Other sizes can be customized