

Carbide Lab Press Mold

Item Number: PMW



Introduction

Form ultra-hard samples with Carbide Lab Press Mold. Made of Japanese high-speed steel, it has a long service life. Custom sizes available.

[Learn More](#)

Model	PMW-A	PMW-B	PMW-C
Material	Carbide YT15		
Indenter hardness	HRC85-HRC90		
Sample size	φ3 /φ4 /φ5 /φ6 /φ8 /φ10 mm	φ12 /φ13 /φ15 /φ18 /φ20 mm	φ22 /φ25 /φ28 /φ30 mm
Cavity depth	30 mm	40 mm	45 mm
External dimensions	φ43×93 mm	φ53×120 mm	φ73×133 mm
Weight(Kg)	0.78	1.8	3.8

Other sizes can be customized