

Pulse Vacuum Lifting Sterilizer

Item Number: KT-Q02



Introduction

The pulse vacuum lifting sterilizer is a state-of-the-art equipment for efficient and precise sterilization. It uses pulsating vacuum technology, customizable cycles, and a user-friendly design for easy operation and safety.

[Learn More](#)

Model	KT-60SV	KT-80SV	KT-100SV	KT-120SV
Effective volume of sterilization chamber	60L ϕ 396×490mm	80L ϕ 396×650mm	100L 504×404×502mm	120L 504×404×590mm
Height after lifting	1310	1590	1475	1565
Work pressure	0.22Mpa			
Operating temperature	134°C			
Thermal uniformity	$\leq \pm 1^\circ\text{C}$			
Timing selection range	0-99min			
Temperature selection range	105-134°C			
Power/supply voltage	Steam generator power: 3KW Drying belt power: 2.4KW/AC220V 50HZ		Steam generator power: 3.6KW Drying belt power: 2.4KW/AC220V 50HZ	
Dimensions (mm)	686×556×913	686×556×1035	840×600×970	840×600×1060
Transport volume (mm)	789×656×1013	786×656×1135	940×610×1000	940×610×1075
Gross weight/net weight	210Kg/180Kg	230Kg/200Kg	300Kg/270Kg	325Kg/295Kg