

## 진동 디스크 밀 소형 실험실 연삭기

품목 번호: KT-JZ



### 소개

효율적인 실험실 연삭을 위한 다용도 진동 디스크 밀에 대해 알아보세요. 지질학, 야금학, 생물학 등에 이상적입니다. 지금 살펴보세요!

### 자세히 알아보기

일련 번호	모델	공급(mm)	배출량(메쉬)	적재량(g)	수량	전압(V)	포트 전압(kW)	크기(LxWxH, mm)	무게(kg)
1	JZGJ-50-1		80-300	50	1	220	0.18	300x300x400	45
2	JZGJ-50-1D		80-300	80-300	1	220	0.18	300x300x400	45
3	JZGJ-80-1		80-300	80	1	220	0.25	380x380x600	65
4	JZGJ-80-1D		80-300	80	1	220	0.25	380x380x600	65
5	JZGJ-80-2		80-300	80x2	2	220	0.25	380x380x600	80
6	JZGJ-80-2D		80-300	80x2	2	220	0.25	380x380x600	80
7	JZGJ-100-1		80-300	100	1	380	1.1	540x520x780	130
8	JZGJ-100-1D		80-300	100	1	380	1.1	540x520x780	130
9	JZGJ-100-2		80-300	100x2	2	380	1.1	540x520x780	160
10	JZGJ-100-2D		80-300	100x2	2	380	1.1	540x520x780	160
11	JZGJ-100-3		80-300	100x3	3	380	1.5	550x600x850	200
12	JZGJ-100-3D		80-300	100x3	3	380	1.5	550x600x850	200
13	JZGJ-100-4		80-300	100x4	4	380	2.2	600x550x850	210
14	JZGJ-100-4D		80-300	100x4	4	380	2.2	600x550x850	210
15	JZGJ-100-6		80-300	100x6	6	380	2.2	700x650x900	250
16	JZGJ-100-6D		80-300	100x6	6	380	2.2	700x650x900	250
17	JZGJ-100-7		80-300	100x7	7	380	2.2	700x650x900	260
18	JZGJ-100-7D		80-300	100x7	7	380	2.2	700x650x900	260
19	JZGJ-200-1		80-300	200	1	380	1.5	540x520x780	180
20	JZGJ-200-1D		80-300	200	1	380	1.5	540x520x780	180
21	JZGJ-200-2		80-300	200x2	2	380	1.5	550x600x850	210
22	JZGJ-200-2D		80-300	200x2	2	380	1.5	550x600x850	210

일련 번호	모델	공급(mm)	배출량(메쉬)	적재량(g)	수량	전압(V)	포트 전압(kW)	크기(LxWxH, mm)	무게(kg)
23	JZGJ-200-3		80-300	200x3	3	380	2.2	700x650x900	220
24	JZGJ-200-3D		80-300	200x3	3	380	2.2	700x650x900	220
25	JZGJ-400-1		80-300	400	1	380	1.5	600x550x850	190
26	JZGJ-400-1D		80-300	400	1	380	1.5	600x550x850	190
27	JZGJ-1000-1D		80-300	1000	1	380	2.2	700x650x900	260
28	JZGJ-100-1DB		80-300	100	1	380	1.1	700x520x800	170
29	JZGJ-100-2DB		80-300	100x2	2	380	1.1	700x520x800	180
30	JZGJ-100-3DB		80-300	100x3	3	380	2.2	1000x580x900	210
31	JZGJ-100-4DB		80-300	100x4	4	380	2.2	1000x580x900	240
32	JZGJ-200-1DB		80-300	200	1	380	1.5	700x520x800	190